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Determination of propofol in human serum by HPLC using a semi-micro
bore column and a column-switching technique

Yoshimi HASHIMOTO*, Ryoko CHIBA**, Yoshimi KOBAYASHI***, Hiroyuki MORIYAMA*,
Yoshiaki MATSUMOTO, Akira TANAKA and Mitsuo MATSUMOTO**

* Separation Center, Tosoh Co., 4560, Kaisei-cho, Shinnanyo-shi, Yamaguchi 746

** Showa College of Pharmaceutical Sciences, 3-3165, Higashitamagawagakuen, Machida-shi, Tokyo 194

*** Scientific Instrument Division, Tosoh Co., 2743-1, Hayakawa, Ayase-shi, Kanagawa 252

(Received 25 November 1996, Accepted 16 January 1997)

An analytical method for the determination of propofol (2,6-diisopropylphenol) in human serum was investigated by high-performance liquid chromatography using a semi-micro bore column and a column-switching technique. The semi-micro bore column [TSKgel ODS-80Ts (2 mm i.d. × 150 mm)] gave a higher peak height (ca. 3.5 folds) compared with that of a conventional column (4.6 mm i.d.). As the result of an examination of some pretreatment columns, higher recoveries of propofol on TSK precolumn BSA-ODS and TSKprecolumn PW were obtained. TSKprecolumn PW was adopted as a pretreatment column, because its less adsorptive ability was considered to be effective in eliminating the hydrophobic compounds in human serum. Under the optimized analytical conditions, the linearity of the calibration curve was obtained over the range of 5~300 ng/ml and the detection limit was 1.5 ng/ml ($S/N=3$). The relative standard deviation against a practical sample was 2.4% ($n=6$). As a consequence, this method was suggested to be utilized for a highly sensitive quantitative analysis of propofol in human serum of an anesthetized patient.

Keywords : HPLC; propofol; column switching; semi-micro bore column; human serum.

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Table 1 Comparison of peak heights and widths between 4.6 mm i.d. and 2 mm i.d. column

	Column inside diameter	
	4.6 mm	2 mm
Propofol		
peak height ^{a)}	131.7	454.9
peak width ^{b)}	6.2	11.8
2-t-Butyl-6-methylphenol		
peak height ^{a)}	85.5	302.4
peak width ^{b)}	6.8	11.5

a) integrator output (unit, mV), b) width at peak half height (unit, s)

Table 2 Characteristics of TSKprecolumn series

	BSA-ODS	PW	SW
Base material	silica gel	polymer gel	silica gel
Hydrophobic layer	ODS(C ₁₈)	-CH ₂ -	-CH ₂ -
Hydrophilic layer	BSA	diol	diol
Particle diameter	20 μm	12 μm	10 μm
pH range	2~7.5	2~12	2~7.5



Fig. 2 Effect of pretreatment time on recovery of propofol

Analytical column: TSKgel ODS-80Ts, 2 mm i.d. × 150 mm; Eluent: CH₃CN-10 mM phosphate buffer (pH 6.8)(65 : 35); Flow rate: 0.2 ml/min; Column temperature: 40°C; Pretreatment column: X, TSKprecolumn BSA-ODS; O, TSKprecolumn PW; Δ, TSKprecolumn SW, 4.6 mm i.d. × 35 mm; Buffer for pretreatment: 50 mM phosphate buffer (pH 6.8); Buffer flow rate: 1.5 ml/min; Injectioe volume: 10 μl

**Fig. 3** Chromatogram of standard samples

Pretreatment column: TSKprecolumn PW, 4.6 mm i.d. \times 35 mm; Pretreatment time: 3 min; Peaks: 1, 2-t-butyl-6-methylphenol; 2, propofol; Other conditions were same as Fig. 2.

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
**Fig. 4** Chromatogram of human serum of a patient anesthetized by propofol

Chromographic conditions were same as Fig. 3.

Table 3 Repeatability test ($n = 6$) of the present method

No.	Propofol/ng ml ⁻¹
1	2608
2	2605
3	2757
4	2605
5	2587
6	2640
Mean	2634
RSD(%)	2.4

Sample: human serum of a patient anesthetized by propofol

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