SEPARATION REPORT

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

Titration curve of TSKgel SP-5PW.

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
Chromatogram of a mixture of ovalbumin(first peak) and trypsin inhibitor (second peak) obtained on DEAE-5PW with 40 min linear gradient elution from 0.02 M tris-HCl buffer of pH 8.0 to 0.02 M tris-HCl buffer of pH 8.0 containing 0.5 M NaCl at a flow rate of 1 ml/min.


$$R(OA, TI) = 2(V_{TI} - V_{OA}) / (W_{TI} + W_{OA})$$

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
Titration curve of DEAE-5PW

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
Sample loading dependence of resolution between ovalbumin and trypsin inhibitor on DEAE-5PW.

Chromatogram of a mixture of ribonuclease A (first peak) and α -chymotrypsinogen (second peak) obtained on SP-5PW with 40 min linear gradient elution from 0.02 M phosphate buffer of pH 7.0 to 0.02 M phosphate buffer of pH 7.0 containing 0.5 M NaCl at a flow rate of 1 ml/min.

$$R(RN, CT) = 2(V_{CT} - V_{RN}) / (W_{CT} + W_{RN})$$

Sample loading dependence of resolution between ribonuclease A and α -chymotrypsinogen A on SP-5PW.

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PROTEIN ADSORPTION CAPACITY OF TSKgel DEAE-5PW

Protein	Adsorption capacity (vg/ml)
Thyroglobulin	35
Ferritin	32
Bovine serum albumin	43

in 0.05 M Tris-HCl buffer (pH 8.3)

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PROTEIN ADSORPTION CAPACITY OF TSKgel SP-5PW

Protein	Adsorption capacity (mg/ml)
γ -Globulin	40
Hemoglobin	43
Lysozyme	54

in 0.02 M phosphate buffer (pH 6.0)

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RECOVERY OF PROTEINS FROM TSKgel DEAE-5PW

Protein	Recovery (%)
Thyroglobulin	98
Ferritin	99
γ -Globulin	100
Bovine serum albumin	102
Hemoglobin	96
Ovalbumin	104
β -Lactoglobulin	103
Trypsin inhibitor	104
Myoglobin	103

Each protein of 0.4 mg was applied to DEAE-5PW column (75 x 7.5 mm I.D.) in 0.02 M tris-HCl buffer (pH 8.5) and the adsorbed protein was desorbed in 0.02 M tris-HCl buffer (pH 8.5) containing 0.5 M NaCl.

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RECOVERY OF PROTEINS FROM TSKgel SP-5PW

Protein	Recovery (%)
γ -Globulin	98
Hemoglobin	96
Trypsinogen	101
α -Chymotrypsinogen A	98
α -Chymotrypsin	104
Myoglobin	88
Lysozyme	95
Ribonuclease A	100
Cytochrome C	103

Each protein of 0.4 mg was applied to SP-5PW column (75 x 7.5 mm I.D.) in 0.02 M phosphate buffer (pH 6.0) and the adsorbed protein was desorbed in 0.02 M phosphate buffer (pH 6.0) containing 0.5 M NaCl.

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RECOVERY OF ENZYMATIC ACTIVITY
FROM TSKgel DEAE-5PW

Enzyme	Recovery (%)
Catalase	80
Lipoxidase	95
α -Chymotrypsin	93

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RECOVERY OF ENZYMATIC ACTIVITY
FROM TSKgel SP-5PW

Enzyme	Recovery(%)
Lipoxidase	84
α -Chymotrypsin	97

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
Change in Ion-Exchange Capacity of DEAE-5PW and SP-5PW during the Exposure in Alkaline or Acid Solution for 10 Days at 25°C.


Ion-exchanger	Solution	Ion-exchange capacity(meq/ml)	
		Before exposure	After exposure
DEAE-5PW	0.5 N NaOH	0.09	0.09
DEAE-5PW	0.5 N HCl	0.09	0.09
DEAE-5PW	20 % CH ₃ COOH	0.09	0.08
SP-5PW	0.5 N NaOH	0.13	0.13
SP-5PW	0.5 N HCl	0.13	0.13
SP-5PW	20 % CH ₃ COOH	0.13	0.13

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
Separation of commercial lipoxidase by high-performance ion-exchange chromatography
Column: DEAE-5PW Sample loading: 1 mg in 0.1 ml
Elution: 60 min linear gradient from 0.02 M ethanolamine-HCl buffer of pH 9.0 to 0.02 M ethanolamine-HCl buffer of pH 9.0 containing 0.5 M NaCl
Flow rate: 1 ml/min Temperature: 25°C Detection: UV absorbance at 280 nm
Recovery of enzymatic activity was 95 %.


Separation of commercial lipoxidase by high-performance ion-exchange chromatography.
Column: SP-5PW Sample loading: 1 mg in 0.1 ml
Elution: 60 min linear gradient from 0.02 M acetate buffer of pH 4.5 to 0.02 M acetate buffer of pH 4.5 containing 0.5 M Na₂SO₄
Flow rate: 1 ml/min Temperature: 25°C Detection: UV absorbance at 280 nm
Recovery of enzymatic activity was 84 %.

Separation of commercial trypsinogen by high-performance ion-exchange chromatography
Column: DEAE-5PW Sample loading: 0.3 mg in 0.1 ml
Elution: 60 min linear gradient from 0.02 M 1,3-diaminopropane-HCl buffer of pH 9.8 to 0.02 M 1,3-diaminopropane-HCl buffer of pH 9.8 containing 0.5 M NaCl
Flow rate: 1 ml/min Temperature: 25°C Detection: UV absorbance at 280 nm

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
Separation of trypsinogen by high-performance ion-exchange chromatography

Column: SP-5PW Sample loading: 0.2 mg
 Elution: 90 min linear gradient from 0.02 M phosphate buffer(pH 7.0)
 to 0.02 M phosphate buffer containing 0.5 M NaCl(pH 7.0)
 Flow rate: 1 ml/min Temperature: 25°C
 Detection: UV absorbance at 280 nm

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
Separation of commercial catalase by high-performance ion-exchange chromatography

Column: DEAE-5PW Sample loading: 1 mg in 0.1 ml
 Elution: 60 min linear gradient from 0.02 M ethanolamine-HCl buffer of pH 9.0
 to 0.02 M ethanolamine-HCl buffer of pH 9.0 containing 0.5 M NaCl
 Flow rate: 1 ml/min Temperature: 25°C Detection: UV absorbance at 280 nm
 Recovery of enzymatic activity was 80 %.

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
Separation of commercial α -chymotrypsin by high-performance ion-exchange chromatography

Column: DEAE-5PW Sample loading: 0.2 mg in 0.1 ml
 Elution: 60 min linear gradient from 0.02 M 1,3-diaminopropane-HCl buffer of pH 9.8
 to 0.02 M 1,3-diaminopropane-HCl buffer of pH 9.8 containing 0.5 M NaCl
 Flow rate: 1 ml/min Temperature: 25°C Detection: UV absorbance at 280 nm
 Recovery of enzymatic activity was 93 %.

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Separation of commercial urease by high-performance ion-exchange chromatography

Column: DEAE-5PW Sample loading: 1 mg in 0.1 ml
 Elution: 60 min linear gradient from 0.02 M tris-HCl buffer of pH 8.0
 to 0.02 M tris-HCl buffer of pH 8.0 containing 0.5 M NaCl
 Flow rate: 1 ml/min Temperature: 25°C Detection: UV absorbance at 280 nm

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