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D. W. Woodall et al., Non-targeted characterization of attributes affecting antibody-FcγRIIIa V158 (CD16a) binding via online affinity chromatography-mass spectrometry, mAb, 14 (2022) No.1, e2004982, CC BY NC  
<https://doi.org/10.1080/19420862.2021.2004982>

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
Figure 1. Overall analytical workflow of the native affinity LC-MS analyses and fraction collection for LC-MS peptide mapping.

Ref.; D. W. Woodall et al., mAb, 14 (2022) No.1, e2004982, CC BY NC

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B. Yao et al., Drifts in N-Linked Glycosylation Result in ADCC Potency Variation of Perjeta® from August 2020 to October 2021 in China, BioMed Research International, 2022, Article ID 7868391, CC BY  
<https://doi.org/10.1155/2022/7868391>

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
FIGURE 2: Drifts in N-linked glycosylation are illustrated, and the reference products can be divided into three groups according to the drifts in %afucose and % galactose. (a) N-glycan profile of the three typical batches of Perjeta®. (b) Scatter plot shows the drifts in %afucose and (c) %galactose. FcγRIIIa affinity analysis further proved drifts in the N-glycan. (d) FcγRIIIa affinity chromatogram of three typical batches of Perjeta®. (e) The percentage of peak I area shows an opposite trend to %galactose. The contents of %afucose and %galactose negatively correlated with the percentage of peak I area, and %galactose exhibited a stronger correlation with FcγRIIIa affinity. (f) Correlation analysis between the percentage of peak I area and the %galactose. (g) Correlation analysis between the percentage of peak I area and the N-glycan (%galactose + %afucose).

Ref.; B. Yao et al., BioMed Research International, 2022, Article ID 7868391, CC BY

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