



# SEPARATION REPORT

## High-performance SEC Semi-Micro Columns Compatible with both Aqueous SEC and Polar Organic SEC, TSKgel SuperAW Series

### Table of Contents

|  |    |
|--|----|
| 1. Introduction                                  | 1  |
| 2. Properties of SuperAW Series                  | 2  |
| 3. Basic Properties of SuperAW Series            | 3  |
| 3-1. Solvent Compatibility                       | 3  |
| 3-2. Calibration Curves with Various Solvents    | 3  |
| 3-3. Resolution and Solvent Consumption          | 5  |
| 3-4. Effect of Flow Rate                         | 6  |
| 3-5. Effect of Sample Injection Volume           | 7  |
| 4. Applications                                  | 8  |
| 4-1. Applications of Various Polymer Measurement | 8  |
| 4-2. Applications in Non-SEC Mode                | 13 |
| 5. Conclusion                                    | 14 |

## 1. Introduction

Size exclusion chromatography (SEC) is a method that separates compounds dissolved in aqueous solution or organic solvents by their molecular sizes. There is a wide range of applicable substances from low molecular weight to high molecular weight and from hydrophilic compounds to hydrophobic compounds, and it is used in various fields. Since it enables not only separation and purification, but also measuring molecular weight or molecular weight distribution of macromolecules, it has especially become an inevitable method in research of macromolecules.

In succession to its organic solvent SEC columns, TSKgel H<sub>XL</sub> series, Tosoh has commercialized several groups of products, including H<sub>HR</sub> series which allows the use of various organic solvents and SuperH column series for ultra-fast GPC. Meanwhile, PW<sub>XL</sub> series has been commercialized as a group of aqueous SEC columns and  $\alpha$  series as a group of SEC columns that can be used widely with aqueous solutions to polar organic solvents.

It has recently commercialized TSKgel SuperAW series, which has been developed by packing hydrophilic base materials with small particle size in semi-micro column and allows solvent replacement from aqueous solutions to organic solvents.

In this article, the basic features of TSKgel SuperAW series and applications on various samples are introduced.

**Table-1 Features of TSKgel Super AW**

| Feature   | Advantage   |
|---|---|
| 1) Micro-particle gel packed in semi-micro column   | <ul style="list-style-type: none"><li>• SEC measurement with short time and high resolution possible<br/>→ Resolution equivalent to PW<sub>XL</sub> or <math>\alpha</math> column (30cm) with nearly half the analysis time</li><li>• Solvent-saving<br/>→ Solvent consumption approximately 1/3 of conventional columns (30cm)</li></ul> |
| 2) Hydrophilic gel with small swelling or shrinkage | <ul style="list-style-type: none"><li>• Solvent replacement from water to organic solvent possible</li><li>• Small sample absorption in polar organic solvent</li><li>• SEC possible in aqueous systems to polar organic solvent systems</li></ul>  |
| 3) High mechanical strength                         | <ul style="list-style-type: none"><li>• High durability</li></ul>   |

## 2. Properties of SuperAW Series

The polymer-type packing materials that have been developed up to present had large swelling or shrinkage properties depending on the difference in solvents used as eluent, and they may have interaction with samples according to its properties when specific packing material is used for various samples with different properties. Therefore, it was necessary that columns were selected for each type of eluent or sample to be analyzed. TSKgel  $\alpha$  series was introduced as a column to solve these problems, which has been received well as a SEC column that covers from aqueous systems to polar organic solvent systems. TSKgel SuperAW series that has been newly developed is capable of achieving separation equivalent to conventional column with 1/2 for analysis time and 1/3 for solvent consumption, due to the micro-particle packing material that has been packed into

semi-micro column. In addition, since hydrophilic base material has been used as the packing material as in  $\alpha$  series, suppression of hydrophobic interaction with sample and SEC measurement under polar organic solvents became possible. With drastic improvements in swelling/shrinkage properties of the packing material, solvent can be replaced in a wide range from water (aqueous solutions) to organic solvents. The features of TSKgel SuperAW series are summarized in Table-1. TSK-Gel SuperAW series consists of 5 types of columns in which gels with different separation ranges and 1 type of mixed column. These 6 types of columns cover a wide range of molecular weights, enabling selection of a column most suited to sample molecular weight or purpose of measurement. Table-2 shows the list of TSKgel SuperAW series products.

**Table-2 List of TSKgel SuperAW Series**


| Grade              | Exclusion limit<br>(PEO/DMF) | Particle size<br>( $\mu\text{m}$ ) | Theoretical plates | Column size<br>(mm I.D. x cm) |
|--------------------|------------------------------|------------------------------------|--------------------|-------------------------------|
| TSKgel SuperAW2500 | $2 \times 10^3$              | 4                                  | > 16,000           | 6.0 x 15                      |
| TSKgel SuperAW3000 | $6 \times 10^4$              | 4                                  | > 16,000           | 6.0 x 15                      |
| TSKgel SuperAW4000 | $4 \times 10^5$              | 6                                  | > 10,000           | 6.0 x 15                      |
| TSKgel SuperAW5000 | $4 \times 10^6$              | 7                                  | > 10,000           | 6.0 x 15                      |
| TSKgel SuperAW6000 | $> 4 \times 10^7$            | 9                                  | > 6,000            | 6.0 x 15                      |
| TSKgel SuperAWM-H  | $> 4 \times 10^7$            | 9                                  | > 6,000            | 6.0 x 15                      |

I.D. indicates the internal diameter.

### 3. Basic Properties of SuperAW Series

#### 3-1. Solvent Compatibility

Figure-1 shows the results of solvent compatibility test on TSKgel SuperAW series. In this test, theoretical plates were measured after substituting from water to test solvent (flow rate: 0.6mL/min, 5 hours or longer), standing (14 hours or longer), and then substituting with water (flow rate: 0.6mL/min) again. This was repeatedly performed on various test solvents. It is evident that in the column of any grade, change in theoretical plates was small even after exchanging with various solvents and that they allow replacement of various solvents from aqueous solutions to polar organic solvents.




**Figure-1 Solvent Compatibility of TSKgel SuperAW Series**

(Theoretical Plates Measurement Conditions)

Column: TSKgel SuperAW Series  
(6.0mm I.D. × 15cm)  
Eluent: Water  
Flow rate: 0.6mL/min  
Temperature: 25°C  
Detection: Refractive index detector  
Sample: Ethylene glycol  
Injection volume: 5μL (2.5g/L)


#### 3-2. Calibration Curves with Various Solvents

Although TSKgel SuperAW series allows exchange to various solvents, it is necessary that the standard sample used should be selected to suit the solvent in actual use. Figures-2 to -6 show calibration curves in water (standard sample PEO, PEG, and pullulan), methanol (PEO, PEG), DMF (PEO, PEG), and DMSO (pullulan). While a calibration curve with favorable linearity has been achieved in each solvent, some difference was seen in molecular weight separation range and slope depending on the solvent type.




**Figure-2 Calibration Curves of TSKgel SuperAW Series (1)**

Column: TSKgel SuperAW Series  
(6.0mm I.D. × 15cm)  
Eluent: Water  
Flow rate: 0.6mL/min  
Temperature: 25°C  
Detection: Refractive index detector  
Samples: Standard polyethylene oxide, polyethylene glycol, ethylene glycol




**Figure-3 Calibration Curves of TSKgel SuperAW Series (2)**

Column: TSKgel SuperAW Series  
(6.0mm I.D. × 15cm)  
Eluent: Water  
Flow rate: 0.6mL/min  
Temperature: 25°C  
Detection: Refractive index detector  
Samples: Standard pullulan, oligomaltose, glucose




**Figure-4 Calibration Curves of TSKgel SuperAW Series (3)**

Column: TSKgel SuperAW Series  
(6.0mm I.D. × 15cm)  
Eluent: Methanol containing 10mmol/L LiBr  
Flow rate: 0.6mL/min  
Temperature: 25°C  
Detection: Refractive index detector  
Samples: Standard polyethylene oxide, polyethylene glycol, ethylene glycol



**Figure-5 Calibration Curves of TSKgel SuperAW Series (4)**

Column: TSKgel SuperAW Series  
(6.0mm I.D. × 15cm)  
Eluent: DMF containing 10mmol/L LiBr  
Flow rate: 0.6mL/min  
Temperature: 25°C  
Detection: Refractive index detector  
Samples: Standard polyethylene oxide, polyethylene glycol, ethylene glycol



**Figure-6 Calibration Curves of TSKgel SuperAW Series (5)**

Column: TSKgel SuperAW Series  
(6.0mm I.D. × 15cm)  
Eluent: DMSO containing 10mmol/L NaNO<sub>3</sub>  
Flow rate: 0.6mL/min  
Temperature: 25°C  
Detection: Refractive index detector  
Samples: Standard pullulan, oligomaltose, glucose


### 3-3. Resolution and Solvent Consumption

Compared to the present TSKgel PW<sub>XL</sub> series or TSKgel  $\alpha$  series, TSKgel SuperAW series has been packed with gel with smaller particle size. The column performance depends on the particle size, and the column efficiency increases as the particle size decreases. Furthermore, it has theoretical plates per unit length that is twice compared to TSKgel PW<sub>XL</sub> series since it is packed in a column with a small dead volume structure. Figure-7 shows chromatograms of dextran T-40 hydrolysate measured on TSKgel SuperAW2500 and TSKgel G2500PW<sub>XL</sub>. It is apparent that SuperAW2500 column

yields separation equivalent to that of the conventional column in approximately half the time.

Moreover, since the flow rate of SuperAW2500 is 60% of G2500PW<sub>XL</sub> column (0.6mL/min compared to 1.0mL/min) in this experiment, the amount of solvent consumption in one sample measurement will be 1/3.

As it is apparent, TSKgel SuperAW series can be called a group of columns with very high cost performance, which is capable of obtaining chromatograms equivalent to those of the conventional general SEC columns with 1/2 measurement time and 1/3 solvent consumption.



**Figure-7 Comparison of Chromatograms**

Column: TSKgel SuperAW2500 (6.0mm I.D. × 15cm)  
TSKgel G2500PWXL (7.8mm I.D. × 30cm)

Eluent: Water

Flow rate: 0.6mL/min (TSKgel SuperAW2500)  
1.0mL/min (TSKgel G2500PWXL)


Temperature: 25°C

Detection: Refractive index detector

Sample: Dextran T-40 hydrolysate

### 3-4. Effect of Flow Rate


Flow rate dependence of HETP (height equivalent to a theoretical plate) was confirmed. Normally, the optimal flow rate depends on the particle size and the molecular weight of the sample. Figure-8 shows the flow rate dependence of HETP on ethylene glycol in SuperAW2500 and G2500PW<sub>XL</sub>. While the optimal linear velocity of PW<sub>XL</sub> column is 10mm/min (flow rate 0.5mL/min) and HETP changes drastically around this value, SuperAW column has smallest HETP in the linear velocity range of 10 to 20mm/min (flow rate 0.3 to 0.6mL/min), yielding almost no change. That is, SuperAW column indicates high theoretical plates in this flow rate



**Figure-8 Relationship between HETP and Linear Velocity (1)**

Column: TSKgel SuperAW2500 (6.0mm I.D. × 15cm)  
 TSKgel G2500PW<sub>XL</sub> (7.8mm I.D. × 30cm)  
 Eluent: Water  
 Temperature: 25°C  
 Sample: Ethylene glycol 2.5g/L  
 Injection volume: 5μL (TSKgel SuperAW2500)  
 10μL (TSKgel G2500PW<sub>XL</sub>)

range and has little deterioration in performance even under high flow rate. Figure-9 shows the HETP flow rate dependence for high-molecular weight pullulan (P-20; molecular weight 23,700, P-5; molecular weight 5,800). In both PW<sub>XL</sub> columns and SuperAW columns, HETP becomes larger as the linear velocity increases. However, it is clear that the rate of change is smaller in SuperAW columns. Based on the above, SuperAW series is evidently a group of columns with extremely small effect of flow rate on sample separation, which can exert high column efficiency in a wide range of flow rates.




**Figure-9 Relationship between HETP and Linear Velocity (2)**

Column: TSKgel SuperAW3000 (6.0mm I.D. × 15cm)  
 TSKgel G3000PW<sub>XL</sub> (7.8mm I.D. × 30cm)  
 Eluent: Water  
 Temperature: 25°C  
 Samples: Pullulan P-5 (1g/L), pullulan P-20 (1g/L)  
 Injection volume: 5μL (TSKgel SuperAW3000)  
 10μL (TSKgel G3000PW<sub>XL</sub>)

### 3-5. Effect of Sample Injection Volume


When the sample injection volume is increased, the peak disperses and affects separation. The effect of injection volume on HETP using a macromolecule (pullulan) with molecular weight distribution and low-molecular weight compound (ethylene glycol) is shown (Figures-10 and -11). While increase in HETP is seen on PW<sub>XL</sub> column when it exceeds 20 $\mu$ L for ethylene glycol, HETP increases from a little beyond 5 $\mu$ L on SuperAW column. On the other hand, no steep increase in HETP is seen

with pullulan until it reaches about 20 $\mu$ L on SuperAW column. Thus it is necessary that sample injection volume of 5 $\mu$ L per column is set when SuperAW column is used for separation of low-molecular weight substances or oligomers, and 20 $\mu$ L per column is set when it is used for measurement of molecular weight or molecular weight distribution of macromolecules.



**Figure-10 Relationship between HETP and Sample Injection Volume (1)**

Column: TSKgel SuperAW2500 (6.0mm I.D.  $\times$  15cm)  
 TSKgel G2500PW<sub>XL</sub> (7.8mm I.D.  $\times$  30cm)  
 Eluent: Water  
 Flow rate: 0.6mL/min (TSKgel SuperAW2500)  
 1.0mL/min (TSKgel G2500PW<sub>XL</sub>)  
 Temperature: 25°C  
 Sample: Ethylene glycol 2.5g/L



**Figure-11 Relationship between HETP and Sample Injection Volume (2)**

Column: TSKgel SuperAW3000 (6.0mm I.D.  $\times$  15cm)  
 TSKgel G3000PW<sub>XL</sub> (7.8mm I.D.  $\times$  30cm)  
 Eluent: Water  
 Flow rate: 0.6mL/min (TSKgel SuperAW3000)  
 1.0mL/min (TSKgel G3000PW<sub>XL</sub>)  
 Temperature: 25°C  
 Sample: Pullulan P-5 (1g/L)  
 Pullulan P-20 (1g/L)  
 Injection volume: 5 $\mu$ L (TSKgel SuperAW3000)  
 10 $\mu$ L (TSKgel G3000PW<sub>XL</sub>)




## 4. Applications

### 4-1. Applications of Various Polymer Measurement

A list of various samples measured on SuperAW column is shown in Table-3. Furthermore, chromatograms of these samples are shown in Figures-12 to -27.


**Table-3 Applications of Various Samples on TSKgel SuperAW Series**

| Figure | Sample name                            | Column            | Solvent                                      |
|--------|--|-------------------|--|
| 12     | Sodium chondroitin sulfate             | TSKgel SuperAWM-H | 0.2mol/L sodium nitrate                      |
| 13     | Sodium alginate                        | TSKgel SuperAWM-H | 0.2mol/L sodium nitrate                      |
| 14     | Carboxymethyl cellulose                | TSKgel SuperAWM-H | 0.2mol/L sodium nitrate                      |
| 15     | Sodium polystyrene sulfonate           | TSKgel SuperAWM-H | 0.2mol/L sodium nitrate/acetonitrile = 80/20 |
| 16     | Polyvinyl pyrrolidone                  | TSKgel SuperAWM-H | 0.2mol/L sodium nitrate/acetonitrile = 80/20 |
| 17     | Gum arabic                             | TSKgel SuperAWM-H | 0.2mol/L sodium nitrate/acetonitrile = 80/20 |
| 18     | Ethylhydroxy-ethylcellulose            | TSKgel SuperAWM-H | Methanol containing 10mmol/L LiBr            |
| 19     | Vinyl alcohol/vinyl butyral copolymer  | TSKgel SuperAWM-H | Methanol containing 10mmol/L LiBr            |
| 20     | Hydroxypropylcellulose                 | TSKgel SuperAWM-H | Methanol containing 10mmol/L LiBr            |
| 21     | Polymethyl vinyl ether                 | TSKgel SuperAWM-H | Methanol containing 10mmol/L LiBr            |
| 22     | Cellulose acetate                      | TSKgel SuperAWM-H | DMF containing 10 mmol/L LiBr                |
| 23     | N-isopropyl acrylamide                 | TSKgel SuperAWM-H | DMF containing 10 mmol/L LiBr                |
| 24     | Polyacrylonitrile                      | TSKgel SuperAWM-H | DMF containing 10 mmol/L LiBr                |
| 25     | Vinyl chloride/vinyl acetate copolymer | TSKgel SuperAWM-H | DMF containing 10 mmol/L LiBr                |
| 26     | Styrene/allyl alcohol copolymer        | TSKgel SuperAWM-H | DMF containing 10 mmol/L LiBr                |
| 27     | Poly (p-phenylene ether sulfone)       | TSKgel SuperAWM-H | DMF containing 10 mmol/L LiBr                |




**Figure-12 Chromatogram of Sodium Chondroitin Sulfate**

Column: TSKgel SuperAWM-H  
(6.0mm I.D. × 15cm × 2)  
Eluent: 0.2mol/L sodium nitrate  
Flow rate: 0.6mL/min  
Temperature: 40°C  
Detection: Refractive index detector  
Sample load: 20μL (0.5g/L)




**Figure-13 Chromatogram of Sodium Alginate**

Column: TSKgel SuperAWM-H  
(6.0mm I.D. × 15cm × 2)  
Eluent: 0.2mol/L sodium nitrate  
Flow rate: 0.6mL/min  
Temperature: 40°C  
Detection: Refractive index detector  
Sample load: 20μL (0.5g/L)




**Figure-14 Chromatogram of Carboxymethyl Cellulose**

Column: TSKgel SuperAWM-H  
(6.0mm I.D. × 15cm × 2)  
Eluent: 0.2mol/L sodium nitrate  
Flow rate: 0.6mL/min  
Temperature: 40°C  
Detection: Refractive index detector  
Sample load: 20μL (0.5g/L)




**Figure-15 Chromatogram of Sodium Polystyrene Sulfonate**

Column: TSKgel SuperAWM-H  
(6.0mm I.D. × 15cm × 2)  
Eluent: 0.2mol/L sodium nitrate/acetonitrile = 80/20  
Flow rate: 0.6mL/min  
Temperature: 40°C  
Detection: Refractive index detector  
Sample load: 20μL (0.5g/L)




**Figure-16 Chromatogram of Polyvinyl Pyrrolidone**

Column: TSKgel SuperAWM-H  
(6.0mm I.D. × 15cm × 2)  
Eluent: 0.2mol/L sodium nitrate/acetonitrile = 80/20  
Flow rate: 0.6mL/min  
Temperature: 40°C  
Detection: Refractive index detector  
Sample load: 20μL (0.5g/L)




**Figure-17 Chromatogram of Gum Arabic**

Column: TSKgel SuperAWM-H  
(6.0mm I.D. × 15cm × 2)  
Eluent: 0.2mol/L sodium nitrate/acetonitrile = 80/20  
Flow rate: 0.6mL/min  
Temperature: 40°C  
Detection: Refractive index detector  
Sample load: 20μL (0.5g/L)




**Figure-18 Chromatogram of Ethylhydroxy-ethylcellulose**

Column: TSKgel SuperAWM-H  
(6.0mm I.D. × 15cm × 2)  
Eluent: Methanol containing 10mmol/L LiBr  
Flow rate: 0.6mL/min  
Temperature: 40°C  
Detection: Refractive index detector  
Sample load: 20µL (0.5g/L)




**Figure-19 Chromatogram of Vinyl Alcohol/vinyl Butyral Copolymer**

Column: TSKgel SuperAWM-H  
(6.0mm I.D. × 15cm × 2)  
Eluent: Methanol containing 10mmol/L LiBr  
Flow rate: 0.6mL/min  
Temperature: 40°C  
Detection: Refractive index detector  
Sample load: 20µL (0.5g/L)




**Figure-20 Chromatogram of Hydroxypropylcellulose**

Column: TSKgel SuperAWM-H  
(6.0mm I.D. × 15cm × 2)  
Eluent: Methanol containing 10mmol/L LiBr  
Flow rate: 0.6mL/min  
Temperature: 40°C  
Detection: Refractive index detector  
Sample load: 20µL (0.5g/L)




**Figure-21 Chromatogram of Polymethyl Vinyl Ether**

Column: TSKgel SuperAWM-H  
(6.0mm I.D. × 15cm × 2)  
Eluent: Methanol containing 10mmol/L LiBr  
Flow rate: 0.6mL/min  
Temperature: 40°C  
Detection: Refractive index detector  
Sample load: 20µL (0.5g/L)




**Figure-22 Chromatogram of Cellulose Acetate**

Column: TSKgel SuperAWM-H  
 (6.0mm I.D. × 15cm × 2)  
 Eluent: DMF containing 10mmol/L LiBr  
 Flow rate: 0.6mL/min  
 Temperature: 40°C  
 Detection: Refractive index detector  
 Sample load: 20μL (0.5g/L)




**Figure-23 Chromatogram of N-isopropyl Acrylamide**

Column: TSKgel SuperAWM-H  
 (6.0mm I.D. × 15cm × 2)  
 Eluent: Methanol containing 10mmol/L LiBr  
 Flow rate: 0.6mL/min  
 Temperature: 40°C  
 Detection: Refractive index detector  
 Sample load: 20μL (0.5g/L)




**Figure-24 Chromatogram of Polyacrylonitrile**

Column: TSKgel SuperAWM-H  
 (6.0mm I.D. × 15cm × 2)  
 Eluent: DMF containing 10mmol/L LiBr  
 Flow rate: 0.6mL/min  
 Temperature: 40°C  
 Detection: Refractive index detector  
 Sample load: 20μL (0.5g/L)




**Figure-25 Chromatogram of Vinyl Chloride/vinyl Acetate Copolymer**

Column: TSKgel SuperAWM-H  
 (6.0mm I.D. × 15cm × 2)  
 Eluent: DMF containing 10mmol/L LiBr  
 Flow rate: 0.6mL/min  
 Temperature: 40°C  
 Detection: Refractive index detector  
 Sample load: 20μL (0.5g/L)



**Figure-26 Chromatogram of Styrene/allylalcohol Copolymer**

Column: TSKgel SuperAWM-H  
 (6.0mm I.D. × 15cm × 2)  
 Eluent: DMF containing 10mmol/L LiBr  
 Flow rate: 0.6mL/min  
 Temperature: 40°C  
 Detection: Refractive index detector  
 Sample load: 20μL (0.5g/L)




**Figure-27 Chromatogram of Poly (p-phenylene Ether Sulfone)**

Column: TSKgel SuperAWM-H  
 (6.0mm I.D. × 15cm × 2)  
 Eluent: DMF containing 10mmol/L LiBr  
 Flow rate: 0.6mL/min  
 Temperature: 40°C  
 Detection: Refractive index detector  
 Sample load: 20μL (0.5g/L)


## 4-2. Applications in Non-SEC Mode

Since it has excellent solvent compatibility, it is possible to obtain different chromatograms as shown in Figure-28 by changing the eluent on samples such as surfactant. That is, the sample is separated based on the molecular size (SEC mode) in 60% acetonitrile solution, and it is retained in other eluent compositions (non-SEC mode). Thus it is possible to set up the elution conditions to suit the purpose of measurement (molecular weight measurement, assay, separation) in one column. Applications of measuring formulated drugs are shown in Figures-29 and -30. In Figure-29, the low-molecular weight components of poultice are retained in the column and clearly separated. In Figure-30, additives are retained in the column and separated. It is apparently very effective in separation of samples containing from low-molecular weight to high-molecular weight components or several low-molecular weight components.




**Figure-28 Measurement of a Surfactant**

Column: TSKgel SuperAW2500  
(6.0mm I.D. × 15cm)  
Eluent: Acetonitrile, acetonitrile solution  
Flow rate: 0.6mL/min  
Temperature: 40°C  
Detection: UV (280nm)  
Injection volume: 20μL



**Figure-29 Application on Poultice**

Column: TSKgel SuperAW2500  
(6.0mm I.D. × 15cm × 2)  
Eluent: Methanol/water = 60/40  
Flow rate: 0.6mL/min  
Temperature: 40°C  
Detection: UV (220nm), refractive index detector  
Injection volume: 10μL



**Figure-30 Application on Cream**

Column: TSKgel SuperAW2500  
 (6.0mm I.D. × 15cm × 2)  
 Eluent: Ethanol  
 Flow rate: 0.6mL/min  
 Temperature: 40°C  
 Detection: UV (275nm)  
 Injection volume: 10μL

## 5. Conclusion

As described in the above sections, it should be clear that TSKgel SuperAW series is a group of columns that covers a wide range from aqueous SEC to polar organic solvent SEC, and has high general versatility achieving high speed, high resolution, and solvent saving compared to conventional columns. In addition, measurement in non-SEC mode is easily examined, and it can be called a first choice column for measuring unknown samples.