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一般名（たんぱく質名および薬物名）	製品名（国際一般名）	会社名	精製プロセス*（推定）	対応可能なTOYOPEARL充填剤
Humanized anti-CD33 IgG4 conjugated to calicheamicin	Mylotarg® (Gemtuzumab ozogamicin)	Pfizer/Wyeth	(AFC-AEC)-SEC	(TOYOPEARL AF-rProtein A HC-650F - GigaCap DEAE-650) - TOYOPEARLHW-55F
Chimeric anti-CD30 IgG1 conjugated to MMAE	Adcetris® (Brentuximab vedotin)	Seattle Genetics/ 武田薬品/Millennium	(AFC-? 3工程)	(TOYOPEARL AF-rProtein A HC-650F)
Humanized anti-HER2 IgG1 conjugated to maytansinoid, DM1	Kadcyla® (Trastuzumab emtansine)	Genentech/Roche	(AFC-CEC)	(TOYOPEARL AF-rProtein A HC-650F - GigaCap S-650)
Humanized anti-CD22 IgG4 conjugated to N-acetyl γ calicheamicin	Besponsa® (Inotuzumab ozogamicin)	Pfizer	(AFC-?)-HIC	(TOYOPEARL AF-rProtein A HC-650F) - Ether-650M, PPG-600M
Humanized anti-CD79b IgG1 conjugated to MMAE	Polivy® (Polatuzumab vedotin-piiq)	Genentech/Roche	(AFC-HIC)-CEC?	(TOYOPEARL AF-rProtein A HC-650F - Phenyl-650) - GigaCap S-650?
Human anti-Nectin 4 IgG1 conjugated to MMAE	Padcev® (Enfortumab vedotin-efyv)	アステラス/ Seattle Genetics	(AFC)	(TOYOPEARL AF-r Protein A HC-650F/AF-Chelate-650M)
Humanized anti-HER2 IgG1 conjugated to DXd	Enhertu® (Trastuzumab deruxtecan-nxky)	AstraZeneca/ 第一三共	(AFC-?)	(TOYOPEARL AF-rProtein A HC-650F - ?)
Humanized anti-TROP2 IgG1 conjugated to SN-38	Trodelyv® (Sacituzumab govitecan-hziy)	Immunomedics	(AFC-?)-HIC	(TOYOPEARL AF-rProtein A HC-650F - ?) - Phenyl-650, PPG-600M
Humanized anti-BCMA IgG1 conjugated to MMAF	Blenrep® (Belantamab mafodotin-blmf)	GlaxoSmithKline	(AFC-AEC)	(TOYOPEARL AF-rProtein A HC-650F - NH2-750F, GigaCap Q-650)
Chimeric anti-CD19 IgG1 conjugated to PBD dimer	Zynlonta® (Loncastuximab tesirine-lpyl)	ADC Therapeutics	(AFC-?) -SEC	(TOYOPEARL AF-rProtein A HC-650F - ?) - HW-55
Human anti-tissue factor IgG1 conjugated to MMAE	Tivdak™ (Tisotumab vedotin-tftv)	Seagen/Genmab	(AFC-?)	(TOYOPEARL AF-rProtein A HC-650F)
Humanized anti-HER2 IgG1 conjugated to MMAE	Aidixi®, 愛地希® (Disitamab vedotin)	RemeGen	(AFC-CEC-SEC)- SEC	(TOYOPEARL AF-rProtein A HC-650F - GigaCap S-650 - HW-55) - HW-40

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1. Y. Matsuda et al., A purification strategy utilizing hydrophobic interaction chromatography to obtain homogeneous species from a site-specific antibody drug conjugate produced by AJICAP™ first generation antibodies, 2020, 9, 16, <https://doi.org/10.3390/antib9020016>
2. S. J. Walsh et al., Site-selective modification strategies in antibody-drug conjugates, Chem. Soc. Rev., 2021, 50, 1305, <https://doi.org/10.1039/D0CS00310G>
3. A. Q. Dean et al., Targeting cancer with antibody-drug conjugates; Promises and challenges, MAbs. 2021; 13(1): 1951427, <https://dx.doi.org/10.1080%2F19420862.2021.1951427>
4. J. T. W. Tong et al., An insight into FDA approved antibody-drug conjugates for cancer therapy, Molecules, 2021, 26, 5847, <https://doi.org/10.3390/molecules26195847>
5. European Medicines Agencies (EMA)の各種医薬品に関する Assessment Report、Scientific Discussion および関連特許

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1. Y. Y. Syed, Tagraxofusp: First global approval, Drugs (2019)79: 579-583, <https://doi.org/10.1007/s40265-019-01087-z>
2. H. Peng, Perspectives on the development of antibody-drug conjugates targeting ROR1 for hematological and solid cancers, Antibody Therapeutics, 2021, Vo. 4 (2021) 222-227, <https://doi.org/10.1093/abt/tbab023>
3. F. Shafiee et al., Targeted Diphtheria Toxin-Based Therapy: A Review Article, Front. Microbiol., 18 October 2019, 1-23 <https://doi.org/10.3389/fmicb.2019.02340>
4. European Medicines Agencies (EMA)の各種医薬品に関する Assessment Report、Scientific Discussion および関連特許

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