

Ordering Information - TSKgel Anion Exchange columns

Part #	Description	Matrix	Housing	ID (mm)	Length (cm)	
19685	TSKgel BioAssist Q, 10 µm, 400 nm	polymer	PEEK	4.6	5	
21410	TSKgel BioAssist Q, 13 µm, 400 nm	polymer	PEEK	10	10	
18761	TSKgel DEAE-2SW, 5 µm 12.5 nm	silica	Stainless Steel	2	25	
07168	TSKgel DEAE-2SW, 5 µm, 12.5 nm	silica	Stainless Steel	4.6	25	
42154	TSKgel Guard Cartridge for 2 mm ID TSKgel DEAE-2SW column, 3 pk, 5 µm	silica	Stainless Steel	2	1	
19308	TSKgel Guard Cartridge Holder for 2 mm ID cartridges		Stainless Steel	2	1	
07648	TSKgel Guardgel Kit for 4.6 mm ID TSKgel DEAE-2SW column, 10 µm	silica	Stainless Steel			
07163	TSKgel DEAE-3SW, 10 µm, 25 nm	silica	Stainless Steel	7.5	7.5	
07648	TSKgel Guardgel Kit for 7.5 mm ID TSKgel DEAE-3SW column, 10 µm	silica	Stainless Steel			
13061	TSKgel DEAE-5PW Glass, 10 µm, 100 nm	polymer	Glass	5	5	
08802	TSKgel DEAE-5PW Glass, 10 µm, 100 nm	polymer	Glass	8	7.5	
14016	TSKgel DEAE-5PW Glass, 13 µm, 100 nm	polymer	Glass	20	15	
18757	TSKgel DEAE-5PW, 10 µm, 100 nm	polymer	Stainless Steel	2	7.5	
07164	TSKgel DEAE-5PW, 10 µm, 100 nm	polymer	Stainless Steel	7.5	7.5	
07574	TSKgel DEAE-5PW, 13 µm, 100 nm	polymer	Stainless Steel	21.5	15	
07930	TSKgel DEAE-5PW, 20 µm, 100 nm	polymer	Stainless Steel	55	20	
08806	TSKgel Glass Guardgel Kit for 5 mm ID and 8 mm ID TSKgel DEAE-5PW columns, 20 µm	polymer	Glass			
14466	TSKgel Glass Guard Column for 20 mm ID TSKgel DEAE-5PW Glass column, 13 µm	polymer	Glass	20	2	
42152	TSKgel Guard Cartridge for 2 mm ID TSKgel DEAE-5PW column, 3 pk, 10 µm	polymer	Stainless Steel	2	1	
19308	TSKgel Guard Cartridge Holder for 2 mm ID cartridges		Stainless Steel	2	1	
07210	TSKgel Guardgel Kit for 7.5 mm ID TSKgel DEAE-5PW column, 20 µm	polymer	Stainless Steel			
16092	TSKgel Guardgel Kit for 21.5 mm ID TSKgel DEAE-5PW column, 20 µm	polymer	Stainless Steel			
07928	TSKgel Guard Column for 55 mm ID TSKgel DEAE-5PW column, 20 µm	polymer	Stainless Steel	45	5	
13075	TSKgel DEAE-NPR, 2.5 µm, nonporous	polymer	Stainless Steel	4.6	3.5	
17088	TSKgel Guard Column for 4.6 mm ID TSKgel DEAE-NPR column, 5 µm	polymer	Stainless Steel	4.6	0.5	
18249	TSKgel DNA-NPR, 2.5 µm, nonporous	polymer	Stainless Steel	4.6	7.5	
18253	TSKgel Guard Column for 4.6 mm ID TSKgel DNA-NPR column, 5 µm	polymer	Stainless Steel	4.6	0.5	



Part #	Description	Matrix	Housing	ID (mm)	Length (cm)	
21962	TSKgel DNA-STAT, 5 µm, nonporous	polymer	Stainless Steel	4.6	10	
07166	TSKgel QAE-2SW, 5 µm, 12.5 nm	silica	Stainless Steel	4.6	25	
07646	TSKgel Guardgel Kit for 4.6 mm ID TSKgel QAE-2SW column, 5 µm	silica	Stainless Steel			
21960	TSKgel Q-STAT, 10 µm, nonporous	polymer	Stainless Steel	3	3.5	
21961	TSKgel Q-STAT, 7 µm, nonporous	polymer	Stainless Steel	4.6	10	
07157	TSKgel SAX, 5 µm, 6 nm	polymer	Stainless Steel	6	15	
08640	TSKgel Sugar AXG, 10 µm, 6 nm	polymer	Stainless Steel	4.6	15	
08639	TSKgel Sugar AXI, 8 µm, 6 nm	polymer	Stainless Steel	4.6	15	
18386	TSKgel SuperQ-5PW Glass, 10 µm, 100 nm	polymer	Glass	8	7.5	
18257	TSKgel SuperQ-5PW, 10 µm, 100 nm	polymer	Stainless Steel	7.5	7.5	
18387	TSKgel SuperQ-5PW, 13 µm, 100 nm	polymer	Stainless Steel	21.5	15	
18388	TSKgel Guardgel Kit for 7.5 mm ID TSKgel SuperQ-5PW column, 20 µm	polymer	Stainless Steel			

Ordering Information - TSKgel Cation Exchange columns

Part #	Description	Matrix	Housing	ID (mm)	Length (cm)	
19686	TSKgel BioAssist S, 7 µm, 130 nm	polymer	PEEK	4.6	5	
21411	TSKgel BioAssist S, 13 µm, 130 nm	polymer	PEEK	10	10	
07167	TSKgel CM-2SW, 5 µm, 12.5 nm	silica	Stainless Steel	4.6	25	
07650	TSKgel Guardgel Kit for TSKgel CM-2SW column, 10 µm	silica	Stainless Steel			
07162	TSKgel CM-3SW, 10 µm, 25 nm	silica	Stainless Steel	7.5	7.5	
07650	TSKgel Guardgel Kit for TSKgel CM-2SW & CM-3SW columns, 10 µm	silica	Stainless Steel			
14012	TSKgel CM-5PW Glass, 13 µm, 100 nm	polymer	Stainless Steel	20	15	
13068	TSKgel CM-5PW, 10 µm, 100 nm	polymer	Stainless Steel	7.5	7.5	
13069	TSKgel Guardgel Kit for 7.5 mm ID TSKgel CM-5PW column, 20 µm	polymer	Stainless Steel			
21965	TSKgel CM-STAT, 10 µm, nonporous	polymer	Stainless Steel	3	3.5	
21966	TSKgel CM-STAT, 7 µm, nonporous	polymer	Stainless Steel	4.6	10	
16653	TSKgel OApak-A, 5 µm, 100 nm	polymer	Stainless Steel	7.8	30	
16654	TSKgel Guard Column for TSKgel OApak-P, 10 µm	polymer	Stainless Steel	6	4	

Part #	Description	Matrix	Housing	ID (mm)	Length (cm)	
07156	TSKgel SCX, Strong Cation Exchange, 5 μ m (Na ⁺), 6 nm	polymer	Stainless Steel	6	15	
07158	TSKgel SCX, Strong Cation Exchange, 5 μ m (H ⁺), 6 nm	polymer	Stainless Steel	7.8	30	
07165	TSKgel SP-2SW, 5 μ m, 12.5 nm	silica	Stainless Steel	4.6	25	
07644	TSKgel Guardgel Kit for 4.6 mm ID TSKgel SP-2SW column, 5 μ m	silica	Stainless Steel			
13062	TSKgel SP-5PW Glass, 10 μ m, 100 nm	polymer	Glass	5	5	
08803	TSKgel SP-5PW Glass, 10 μ m, 100 nm	polymer	Glass	8	7.5	
14017	TSKgel SP-5PW Glass, 13 μ m, 100 nm	polymer	Glass	20	15	
18758	TSKgel SP-5PW, 10 μ m, 100 nm	polymer	Stainless Steel	2	7.5	
07161	TSKgel SP-5PW, 10 μ m, 100 nm	polymer	Stainless Steel	7.5	7.5	
07575	TSKgel SP-5PW, 13 μ m, 100 nm	polymer	Stainless Steel	21.5	15	
07934	TSKgel SP-5PW, 20 μ m, 100 nm	polymer	Stainless Steel	55	20	
08807	TSKgel Glass Guardgel Kit for 5 mm ID & 8 mm ID TSKgel SP-5PW columns, 20 μ m	polymer	Glass			
42153	TSKgel Guard Cartridge for 2 mm ID TSKgel SP-5PW column, 3 pk, 10 μ m	polymer	Stainless Steel	2	1	
19308	TSKgel Guard Cartridge Holder for 2 mm ID cartridges		Stainless Steel	2	1	
07211	TSKgel Guardgel Kit for 7.5 mm ID TSKgel SP-5PW column, 20 μ m	polymer	Stainless Steel			
16093	TSKgel Guardgel Kit for 21.5 mm ID TSKgel SP-5PW column, 20 μ m	polymer	Stainless Steel			
07932	TSKgel Guard Column for 55 mm ID TSKgel SP-5PW column, 20 μ m	polymer	Stainless Steel	45	5	
13076	TSKgel SP-NPR, 2.5 μ m, nonporous	polymer	Stainless Steel	4.6	3.5	
21963	TSKgel SP-STAT, 10 μ m, nonporous	polymer	Stainless Steel	3	3.5	
21964	TSKgel SP-STAT, 7 μ m, nonporous	polymer	Stainless Steel	4.6	10	